


**ETP**

**...THE FAMILY IS GROWING**



 **INTELSINT**




**Fully enclosed, "state of the art" featured, safe and highly reliable, the Intelsint ETP is the right answer for the day-by-day tissue processing needs of every Histology laboratory.**

**With a processing capacity of 360 specimens, three paraffin containers and an easy to operate HMI (Human Machine Interface) it sets itself at the top of the category.**

**Designed with in mind concepts of specimens safety and operator ergonomics, the ETP is capable to perform also in not favorable environments.**

**Its fully automated RFD (Remote Fill and Drain) device allows the user a safe, clean and easy reagents renewal.**

**Simple and economic service and maintenance is another of its key Features.**



Graphical  
Easy and Intuitive  
HMI  
(Human  
Machine  
Interface)

2 Ports for easy  
and  
automatic  
Reagent Remote  
Fill & Drain  
  
1 Port for  
External Wax  
Dumping

4,5 liters  
reagent tanks

The ETP Sample Processing Chamber  
  
Completely made in high quality Stainless Steel  
Safe and easy to operate lid lock mechanism  
Special, easy to clean, surface treatment.  
Gas piston to keep the lid open.  
Lid open/close sensor.

USB Port  
for  
Data Backup  
and  
Software Automatic  
Upgrades

## ETP - Features Check List

- ✓ Overnight Processing
- ✓ Reagent Management System
- ✓ Remote Fill & Drain
- ✓ 360 samples capacity
- ✓ 15" Touch Screen Operation Control
- ✓ 3 Level Sensors
- ✓ Reagent Partial Load Capacity
- ✓ Wax Cleaning Cycle
- ✓ Data Backup on USB port
- ✓ 18 Process Programs
- ✓ Revers Processing Special Program
- ✓ Process End Time Management
- ✓ Fumes Charcoal Filter
- ✓ Anti-blackout and Power Failures Devices
- ✓ Remote Alarm
- ✓ Remote Monitoring and Servicing
- ✓ Factory Pre-filled Reagents

# ETP



## Process Editing

Simple  
Clear  
Easy  
Intuitive

18 Process available

Drag & Copy Function  
Copy Process Function  
Controlled Field Value

## Process Start and Run

Most Used Process Function  
Easy End of Process Setting

Graphic Process Controls  
Multi-language  
Reagent Level Check  
Partial Reagent Load  
Samples Baskets Check

Info Function to check 'at a glance'  
the most important instrument  
parameters

Reagent Quality Check



## Remote Fill & Drain

To change the reagents  
without removing tanks

Graphic RFD Control  
Reagent Quality Check  
Easy and Intuitive HMI  
RFD connected with RMS

Fast and Fully Automated

2 Ports for Reagents  
1 Port for Paraffin

# INTELSINT ETP - Technical Specification

Sample Processing Chambers	1
Sample baskets	3
Sample Processing Capacity	360
Reagent Tanks - 4,5 Liters capacity	12
Paraffin Containers - 4,5 Liters capacity	3
SPC Level Sensors	3
Reagent Partial Load Capacity	Yes
Charcoal Filters	1
Remote Fill & Drain Ports	2
Wax Dumping Port	1
Color Touch Screen	15"
Processing Programs	18
Purge Reagents	2
Reagents Temperature Range (°C)	20-45
Paraffin Temperature Range (°C)	55-65
Process Start delay Time	Up to 2 weeks
RMS - Reagent Management System	Yes
RFD - Remote Fill & Drain	Yes
EWD - External Wax Dumping	Yes
USB Ports	4
Printer Ports	1
Remote Alarm Port	1
Network Ports (Ethernet)	2
Approvals	CE/IVD/UL
Rating (50/60 Hz)	100-240V 50-60Hz
Available software languages	Chinese, English, German, Italian, French, Spanish, Polish, Czech, Portuguese


**INTELSINT**

INTELSINT srl Via Rivoli 122 - 10090 Villarbasse (TO) - Italy  
 Phone +39 011 9528015 - Fax +39 011 952723  
 info@intelsint.com - www.intelsint.it